 **Johnson & Johnson INTERNATIONAL**
c/o European Logistics Centre
Leonardo Da Vincielaan, 15
BE-1831 Diegem
Belgium

1-877-ETHICON
+1-513-337-6928

PROCEDED™ VENTRAL PATCH

8560905

09/2019





LAB0012529v7

© Ethicon, Inc. 2016

CE 2797

- ar** جهاز شبكة قابلة للامتصاص جزئياً
- cs** Částečně vstřebatelný zdravotnický prostředek ze síťoviny
- da** Delvist resorberbar mesh
- de** Teilresorbierbares Netzimplantat
- el** Μερικώς απορροφήσιμη διάταξη πλέγματος
- en** Partially Absorbable Mesh Device
- es** Dispositivo de malla parcialmente absorbible
- fi** Osittain resorboituva verkko
- fr** Dispositif en treillis partiellement résorbable
- hu** Részben felszívódó háló eszköz
- it** Dispositivo a rete parzialmente assorbibile
- kk** Жартылай сіңірлелін тор құрылғысы
- ko** 부분 흡수성 메쉬 장치
- mk** Делумно ресорптивен мрежен инструмент
- nl** Gedeeltelijk resorbeerbaar mesh-hulpmiddel
- no** Delvis resorberbart nett
- pl** Wyrób z częściowo wchłanialnej siatki
- pt** Dispositivo de rede parcialmente absorvível
- ru** Частично рассасывающееся сетчатое устройство

- sk** Čiastočne vstrebateľné sietkové zariadenie
- sr** Delimično resorptivno medicinsko sredstvo sa mrežicom
- sv** Delvis resorberbar nätenhet
- tr** Kismen Emilebilir Mesh Cihazı
- uk** Виріб із частково розсмоктуваної сітки

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LAB-0012529 | Rev:7
Released: 14 Nov 2019
CO: 100750560
Release Level: 4: Production

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INTRODUCTION

This package insert is designed to provide instructions for use of PROCEED™ Ventral Patch. It is not a comprehensive reference to surgical technique for repair of abdominal wall hernias.

PROCEED™ Ventral Patch is intended for use only by physicians who are trained in the surgical procedures and techniques required for hernia repairs and the implantation of synthetic meshes. The selection of mesh for any given patient is a function of numerous factors including, but not limited to, the patient's past medical and surgical history, current medical condition (i.e., comorbidities), surgical technique, and size and location of the hernia. The physician is advised to consult the medical literature regarding techniques, complications, and adverse reactions before selecting a mesh.

DESCRIPTION

PROCEED™ Ventral Patch is a device with absorbable and nonabsorbable components. The visceral side of the patch including the anchoring straps (see Figure 1) is comprised of an absorbable oxidized regenerated cellulose (ORC) fabric attached to a macroporous nonabsorbable mesh knitted from polypropylene filaments. The polymer of the undyed and dyed (phthalocyanine blue, Color Index No. 74160) polypropylene fiber is identical to the material used for dyed/undyed PROLENE™ Suture. The polypropylene mesh is attached to two absorbable polydioxanone films. The polydioxanone polymer is identical to the polymer used in PDS™ II (polydioxanone) Suture. The parietal side of the patch contains a reinforcement film and a positioning ring composed of polydioxanone. A mesh composed of absorbable polyglactin 910 filaments is attached on top of the polydioxanone polymer reinforcement film. The anchoring straps of the patch are designed to facilitate placement and fixation of the device. A nonabsorbable polyester suture is attached to each anchoring strap, enabling tension to be placed on the patch prior to fixation. The polyester polymer is identical to the material used in ETHIBOND™ Suture. The absorbable ORC layer on the visceral side of the mesh is designed to temporarily separate tissue/organ surfaces during healing thereby reducing the risk of tissue attachment to the mesh.

INDICATIONS

PROCEED™ Ventral Patch is indicated for the repair of abdominal wall hernias and abdominal wall deficiencies that require the addition of a reinforcing material to obtain the desired surgical result.

CONTRAINDICATIONS

None known.

WARNINGS

- Failure to properly follow instructions may result in improper functioning of the device and could lead to injury. Please read all information carefully.
- If this device is used in patients with the potential for growth or tissue expansion (such as infants, children, or women who may become pregnant), the surgeon should be aware that the device will not stretch significantly as the patient grows.
- This device is indicated for abdominal wall hernia repair and not for gynecologic procedures. Gynecologic procedures should be performed using devices indicated for gynecologic repairs.
- It is recommended that the device not be used in a contaminated field, as contamination of the device may lead to infection that may require removal of the device.
- As with any implant, an acute and permanent foreign body response will occur. In some patients, this response can result in one or more of the adverse reactions listed below.
- This device is a permanent implant that is designed to integrate into the tissue. In cases in which the device needs to be removed, in part or in whole, significant dissection may be required.

- Care should be taken not to damage the absorbable ORC side of the device during handling as damage may increase the risk of postsurgical complications, including adhesions. Ensure that the device is intact prior to implantation.
- The ORC side must be oriented toward the viscera of the body. For further information see the Application/Instructions for Use section.
- Insufficient overlap on any side of the defect may increase the risk of postoperative complications including recurrence. Consult Application/Instructions for Use section.
- Insufficient or improper fixation may increase the risk of postoperative complications, including recurrence. Consult Application/Instructions for Use section.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and anyone coming in contact with the device.
- Inspect the device carefully before implantation. Do not use the device if it is damaged.

PRECAUTIONS

- Do not soak or stretch the PROCEED™ Ventral Patch.
- The safety and effectiveness of pre-treating PROCEED™ Ventral Patch with solutions other than saline (e.g., peritoneal instillates, and/or medications) prior to implantation have not been studied.

ADVERSE REACTIONS

Potential adverse reactions are those typically associated with surgically implantable materials, including infection, inflammation, seroma formation, acute or chronic pain, foreign body sensation, hematoma, nerve damage, soft tissue injury, adhesion formation, fistula formation, extrusion/erosion, excessive contraction or shrinkage of the tissue surrounding the mesh, and mesh failure/hernia recurrence.

One or more revision surgeries may be necessary to treat the above-mentioned adverse reactions. Revision surgery may not resolve the adverse reactions and may pose a risk of additional adverse reactions.

APPLICATION / INSTRUCTIONS FOR USE

Correct surface orientation is critical for PROCEED™ Ventral Patch to function as intended. Only the visceral side of the mesh is laminated with an absorbable ORC layer. Therefore, the PROCEED™ Ventral Patch must be placed with ORC layer facing the viscera. Active bleeding should be controlled prior to placement of PROCEED™ Ventral Patch. Only the fixation straps can be trimmed.

Sizing and Placement

PROCEED™ Ventral Patch should not be cut or trimmed except for the fixation straps. In ventral or incisional hernia repair, the mesh should extend at least 3 to 5 cm beyond the margins of the hernia defect, unless the surgeon decides that additional overlap onto healthy tissue is needed.

Procedure steps

Just prior to insertion, dip mesh patch in saline for ease of use and to avoid attachment of tissue during insertion. Fold the mesh patch into a semi-circle with the ORC side facing outward (straps folded in) for insertion into the defect (see Figure 2).

During insertion, be certain to avoid kinking the patch (e.g., by using clamps or extensive folding/rolling etc.). Pull the device with sufficient tension to ensure tight seating against the abdominal wall as demonstrated in Figure 3. If it is necessary to remove the device, it should not be reinserted.

Revision

It is important to ensure that no tissue is trapped between the device and the abdominal wall. It is also important to have a small defect at the anterior fascia through the mesh straps (Figure 4). It is recommended that the end of the strap is cut off at least 6.5 mm (1/4 inch) from the fixation point. It is recommended that nonabsorbable sutures be used to fixate the patch.

PERFORMANCE / ACTIONS

An animal study showed that implantation of PROCEED™ Ventral Patch elicits an inflammatory reaction that stimulates the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The absorbable components are essentially absorbed within 238 days, whereas the polypropylene material is not absorbed. Minimal visceral tissue attachment has been demonstrated in an animal study.

STERILITY

The device is sterilized by irradiation. Do not resterilize. Do not use if the package is opened or damaged.

STORAGE

No special storage conditions required. Do not use after expiry date.

HOW SUPPLIED

PROCEED™ Ventral Patch is available in single-use, sterile packets in a variety of sizes.

TRACEABILITY

A traceability label that identifies the type, size, expiry date and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted.

SYMBOLS USED ON LABELING

 Caution	 Manufacturer	 STERILE  Sterilized using irradiation
 Do not reuse	 Quantity	 Rx only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed healthcare practitioner.
 Do not resterilize	 Catalogue number	
 Do not use if package is damaged.	 Batch code	 CE 2797 CE-mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive 93/42/EEC.
	 Use-by date	

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LAB-0012529 | Rev:7
Released: 14 Nov 2019
CO: 100750560
Release Level: 4. Production

THIS PAGE IS LEFT BLANK INTENTIONALLY.

PPE Specification
Labeling Specification
8560905 Proceed Ventral Patch CE Marked IFU

LAB-0012529 | Rev:7
Released: 14 Nov 2019
CO: 100750560
Release Level: 4. Production

THIS PAGE IS LEFT BLANK INTENTIONALLY.

PPE Specification
Labeling Specification
8560905 Proceed Ventral Patch CE Marked IFU

LAB-0012529 | Rev:7
Released: 14 Nov 2019
CO: 100750560
Release Level: 4. Production

THIS PAGE IS LEFT BLANK INTENTIONALLY.

IFU PRINTING SPECIFICATION SHEET

COVER (FLAT)

19.528"
(496 mm)

6.5"
(165.1 mm)

SADDLE STITCH PAGES HERE

TEXT PAGES

6.5"
(165.1 mm)

5"
(127 mm)

COVER (FOLDED)

6.5"
(165.1 mm)

5"
(127 mm)

Binding Method

Saddle stitch

TITLE	DESCRIPTION		LAB NUMBER	SPECIAL INSTRUCTIONS/COMMENTS	BINDING	COLORS			
PROCEED™ Ventral Patch	Multilingual IFU Booklet		LAB0012529v7	n/a	Saddle stitch	Black			
FLAT SIZE 19.528" x 6.5" 496 mm x 165.1 mm	FOLDED SIZE 5" x 6.5" 127 mm x 165.1 mm	RMC NUMBER 8560905	PAGE COUNT 8 - Cover Pgs 80 - Text Pgs	LANGUAGES ar, cs, da, de, el, en, es, fi, fr, hu, it, kk, ko, mk, nl, no, pl, pt, ru, sk, sr, sv, tr, uk	SELF COVER <input type="checkbox"/>	PLUS COVER <input checked="" type="checkbox"/>	SEALING METHOD n/a	WAFER SEAL <input type="checkbox"/>	
BLEED SIZE .5" (12.7 mm) .125" (3.175 mm) <input type="checkbox"/> <input checked="" type="checkbox"/>	NONE <input type="checkbox"/>	BLEED ALL SIDES <input checked="" type="checkbox"/>	BLEED TOP <input type="checkbox"/>	BLEED RIGHT <input type="checkbox"/>	BLEED LEFT <input type="checkbox"/>	BLEED BOTTOM <input type="checkbox"/>	DRAWING IS NOT TO SCALE: DRAWINGS REFLECT INFORMATION FOR PRODUCTION OF PRINTED PIECES AND DO NOT CONTAIN ACTUAL ARTWORK. This document or data herein or herewith is not to be reproduced, used or disclosed in whole or part without the permission of Ethicon, Inc.		
STOCK Refer to the Material or Finished Goods Specification	ETHICON								