CAPSURE™

Permanent Fixation System

Instructions for Use

Mode d'emploi

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
BAIRD
DAVOL INC.

PRODUCT DESCRIPTION

The CAPSURE™ Permanent Fixation System is a sterile single use device that delivers either 15 or 30 permanent polyetheretherketone (PEEK) and 316L stainless steel fasteners. The shaft on the CAPSURE™ Fixation System is 37 cm in length. The fastener, including the PEEK cap and 316L stainless steel coil, is 4.2 mm in length. The fixation instrument shafts have an outer diameter of 5 mm, and may be used in open procedures or with most 5 mm trocars in laparoscopic procedures. The device includes a fastener gauge located on the top of the hand piece. The gauge will move right to left as the fasteners are deployed that indicates the approximate level of fasteners remaining in the device.

15 count fastener level indicator
starts at midpoint and moves
towards left as fasteners deplete.

30 count fastener level indicator
starts on the right and moves
towards left as fasteners deplete.

- REF 0113215 CAPSURE™ Permanent Fixation System, Straight 37 cm length, 15 Fasteners
- REF 0113230 CAPSURE™ Permanent Fixation System, Straight 37 cm length, 30 Fasteners

INDICATIONS

The CAPSURE™ Permanent Fixation System is indicated for the approximation of soft tissue and fixation of surgical mesh to tissues during open or laparoscopic surgical procedures, such as hernia repair.

CONTRAINDICATIONS

1. This device is not intended for use except as indicated.
2. Do not use this device where hemostasis cannot be verified visually after application.
3. Contraindications associated with laparoscopic and open surgical procedures relative to mesh fixation apply, including but not limited to: fixation of vascular or neural structures and in ischemic or necrotic tissue.
4. Carefully inspect the area in the vicinity of the tissue being fastened to avoid inadvertent penetration of underlying structures such as bone, nerves, vessels, and viscera. Use of the CAPSURE™ Permanent Fixation System in the close vicinity of such underlying structures is contraindicated. For reference, the length of the fastener below the fastener head is 3.2 mm, the fastener head is another 1 mm (total 4.2 mm).
5. This device should not be used in tissues that have a direct anatomic relationship to major vascular or nerve structures. This includes the deployment of tacks in the diaphragm in the vicinity of the pericardium, aorta, or inferior vena cava during diaphragmatic hernia repair.

WARNINGS

1. The CAPSURE™ Permanent Fixation System is intended for Single Use Only - **DO NOT RESTERILIZE**. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient. Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient or end user.
2. Do not use beyond the expiration date on the package.
3. Prior to use, carefully examine package and product to verify neither is damaged and that all seals are intact. Do not use if the package is damaged or open.
4. Verify mechanical and electrical compatibility of devices from different manufacturers prior to using them together in a procedure.
5. Prosthetic should be evaluated for compatibility prior to use.
6. Users should be familiar with surgical procedures and techniques involving permanent materials before employing CAPSURE™ Permanent Fixation System fasteners for wound closure, as the risk of wound dehiscence may vary with the site of application and the material used.
7. As with any implant material, the presence of bacterial contamination may enhance bacterial infectivity. Accepted surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

 After use, the CAPSURE™ Permanent Fixation System may be a potential biohazard. Handle and dispose of in accordance with any local and federal laws regarding medical waste.

en

PRECAUTIONS

1. Please read all instructions before using the CAPSURE™ Permanent Fixation System.
2. Only persons having adequate medical training and familiarity with surgical techniques should perform surgical procedures. Consult the medical literature relative to technique, complications and hazards prior to any surgical procedure.
3. The CAPSURE™ Permanent Fixation System can be used with most 5 mm trocars. Ensure compatibility by inserting the device into the trocar prior to introduction into the patient. The CAPSURE™ Permanent Fixation System should enter and exit the trocar easily without excessive force. The use of too much force could damage the instrument.
4. Adequate counter pressure should be applied on the target area. Avoid placing hand or finger directly over the area where fastener is being deployed to prevent injury.
5. Use caution when applying the CAPSURE™ fastener over or in proximity to underlying bone, vessels, nerves, or viscera. The intended fixation site should be assessed to ensure that while the tissue is compressed the total distance from the surface of the tissue to any underlying structures is greater than the length of the CAPSURE™ fastener.
6. Care should be taken not to use excessive counter pressure as this may damage the tissue, the material being fixated, and/or the device.
7. If the device locks and cannot be separated from a fastener that has been deployed into mesh and/or tissue, rotate the device counterclockwise to free the fastener from the tissue and/or to free the device. If the fastener does not deploy properly, remove the device from the patient and test the device in gauze to ensure proper fastener deployment, otherwise discard the device appropriately and use a new CAPSURE™ Permanent Fixation System. Once proper fastener deployment is confirmed, the device may be reinserted into the patient.

ADVERSE REACTIONS

Adverse reactions and potential complications associated with fixation devices such as the CAPSURE™ Permanent Fixation System may include, but are not limited to the following: hemorrhage, pain, edema and erythema at wound site; septicemia/infection; hernia recurrence/wound dehiscence, erosion and allergic response in patients with known sensitivities to PEEK and metals contained in 316L stainless steel, including chromium, nickel, copper, and iron.

INSTRUCTIONS FOR USE

1. Take the CAPSURE™ Permanent Fixation System out of the sterile package using sterile technique.
2. Bring the prosthesis (mesh) or tissue into position. For tissue approximation, be sure that adequate overlap of tissue exists.
3. Place the tip of the CAPSURE™ Permanent System at the desired location and apply adequate counter pressure. Different types of mesh may require different amounts of counter pressure. Adjust angle and counter pressure appropriately.
4. Compress hand piece trigger in a single, complete and uninterrupted stroke to drive a permanent fastener through the mesh into the tissue. Keep consistent counter pressure on the tip of the device through the entire stroke. Release the trigger allowing it to return completely to its resting position. Repeat this procedure until all required fasteners are deployed.
5. The fasteners should be placed entirely in tissue and the head of the fastener should be flush against the mesh or tissue in order to achieve the best fixation performance. If the fastener head is not flush in mesh or tissue, it is recommended to use a Maryland grasper, or similaratraumatic surgical grasping device, and place the tip of one jaw on the inside of the fastener head and the tip of the other jaw on the outer diameter of the head. To fully seat the fastener, turn the grasper clockwise. If the fastener still does not fully seat, use a Maryland grasper in the same manner to remove the fastener by turning the grasper counter clockwise. Place another fastener in the same vicinity as the removed fastener. If the fastener head becomes dislodged from the fastener coil, remove the coil by grasping the proximal end of the coil and rotating the grasper counterclockwise. Loose fasteners should be retrieved from the abdominal or other cavities if possible. Dispose of loose fasteners in accordance with any local and federal laws regarding sharps disposal requirements.
6. Care should be taken to ensure that the prosthesis (mesh) is adequately fixated to the abdominal wall. If necessary, additional fasteners and/or sutures should be used.
7. After successful deployment of all required fasteners, handle and dispose of in accordance with any local and federal laws regarding medical waste.

STORAGE

Store the CAPSURE™ Permanent Fixation System at room temperature. Avoid prolonged exposure to elevated temperatures. The device is for single use only. Do not use if package is damaged or open.

MRI Information

MR Conditional

The CAPSURE™ Fastener is MR Conditional.

Non-clinical testing demonstrated that the CAPSURE™ Fastener is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 1,000-Gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR system

MRI-Related Heating

In non clinical testing, the CAPSURE™ fastener produced the following temperature rise during MRI performed for 15-min of scanning (i.e. per pulse sequence) in a 3-Tesla MR system (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Healthcare, Milwaukee, WI).

MR system reported, whole body average SAR	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.7-W/kg
Highest temperature change	1.8°C
Temperature scaled to whole body averaged SAR of 4-W/kg	2.5°C

Artifact Information

The maximum artifact size as seen on the gradient echo pulse sequence at 3-Tesla extends approximately 10-mm relative to the size of the shape of the CAPSURE™ fastener.

Bard, Davol and CapSure are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
Copyright © 2014 C. R. Bard, Inc. All Rights Reserved.

	Contents
	Do not use if package is damaged
	15 Permanent Fasteners
	30 Permanent Fasteners
	Do Not Resterilize
	U.S. federal law restricts this device to sale by or on the order of a physician

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Bard, Davol et CapSure sont des marques de commerce et/ou des marques déposées de C. R. Bard, Inc. ou d'une filiale. Copyright © 2014 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2014 C. R. Bard, Inc. Tutti i diritti riservati.
Bard, Davol y CapSure son marcas comerciales y/o marcas comerciales registradas de C. R. Bard, Inc. o de una filial. Copyright © 2014 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2014 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2014 C. R. Bard, Inc. Alle rettigheder forbeholdes.
Copyright © 2014 C. R. Bard, Inc. Kaikki oikeudet pidätetään.