## Exempel på ablationszon

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
**Emprint™**

- [REF] CA15L1 **Short Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone  
15 cm
- [REF] CA20L1 **Standard Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone  
20 cm
- [REF] CA30L1 **Long Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone  
30 cm
- [REF] CA15L2 **Short Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone - Reinforced Antenna  
15 cm
- [REF] CA20L2 **Standard Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone - Reinforced Antenna  
20 cm
- [REF] CA30L2 **Long Percutaneous Antenna with Thermosphere™ Technology**  
Large Spherical Ablation Zone - Reinforced Antenna  
30 cm

- ① Radiating shaft (green)  
② Adjustable depth guide (comes in a separate pouch in the package)  
③ Antenna shaft (1.3G)  
④ Antenna handle  
⑤ Antenna cable connector  
⑥ Clips  
⑦ S Clip  
⑧ Saline bar spike connector  
⑨ Orange tubing that goes into the pump  
⑩ Antenna cable

**Indication for Use**

The Emprint™ Ablation System is intended for use in percutaneous, laparoscopic, and intraoperative coagulation (ablation) of soft tissue, including partial or complete ablation of non-resectable liver tumors.

The Emprint™ Ablation System is not intended for use in cardiac procedures.

**Device Description**

For use with Emprint™ Ablation System with Thermosphere™ Technology.

Not made with natural rubber latex

Do not use if package is opened or damaged.

Type BF applied part

Cable connect/disconnect


**Warning**

Do not resterilize. The product cannot be adequately cleaned and/or sterilized by the user to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize the product may result in infection or product failure risks to the patient.

Do not activate microwave energy until system setup is complete and the radiating section of the antenna is fully inserted into tissue.

When treating shallow lesions, see examples of ablation zones beginning on page 83 to ensure proper insertion depth of the antenna to avoid unintended thermal injury. Consider using an open procedure when the ablation zone may not be completely contained in target tissue.

When performing ablations, ensure the green, radiating section of the antenna is fully inserted into the tissue. A shallower insertion depth than 2.8 cm may lead to unintended thermal injury.

**Warning**

Use caution when creating ablation zones near large vasculature, including, but not limited to, the hepatic artery and portal vein. Similar caution should be exercised near ductal structures, the diaphragm, the colon, and bronchi. Consider other treatment options for lesions located in these areas.

The radiating section of the antenna becomes hot when microwave power is applied. Care should be taken to ensure that the radiating section does not come in contact with unintended tissue, as thermal injury may occur. Users should not touch the antenna radiating section at any time during the application of power. Do not activate microwave power to an exposed antenna.

The antenna cable may become hot during a procedure. Take care when handling the cable. Do not allow prolonged contact with the patient's skin during a procedure as thermal injury could occur.

Do not activate the antenna while it is in contact with metal objects or instruments as unintended thermal injury could occur.

Microwave ablation procedures are not recommended for patients with cardiac pacemakers or other implanted electronic devices. Potential risks have not been evaluated.

Do not defibrillate a patient with an antenna inserted. Completely remove the antenna from the patient before defibrillation. Potential risks have not been evaluated.

Microwave procedures are not recommended for pregnant patients. Potential risks to the patient and/or fetus have not been established.

Smoke generated during electrosurgical procedures can be potentially harmful to patients and the surgical team. Adequately vent the smoke by using a surgical-smoke evacuator or other means.

Do not look into or place the radiating section near eyes or testes when energized as this may result in burns or other injuries to the patient or operator.

**Fire/Explosion Hazard: Do not use the Emprint Ablation System in potentially flammable or explosive environments.**

The following substances will contribute to increased fire and explosion hazards in the operating room:

- Flammable substances (such as flammable anesthetics, alcohol-based skin prepping agents and tinctures)
- Naturally occurring flammable gases that may accumulate in body cavities such as the bowel
- Oxygen enriched atmospheres
- Oxidizing agents (such as nitrous oxide [N<sub>2</sub>O] atmospheres)

Heating associated with microwave energy can provide an ignition source. Observe fire precautions at all times. When using microwave energy in the same room with any of these substances or gases, prevent their accumulation or pooling under surgical drapes, or within the area where microwave procedures are performed.

### Caution

The antennas should only be used by clinicians and staff properly trained in the use of this technology and its associated warnings and cautions. Clinicians should also be skilled/trained in the use of image guidance technology such as ultrasound or computed tomography for the placement of antennas.

Emprint™ antennas are shipped sterile and are to be discarded after a single procedure. Do not attempt to resterilize.

### Storage

- Store antennas away from moisture and direct heat.
- Discard the product if it is damaged, or if sterile packaging is damaged or opened.
- Discard products that exceed the expiration date displayed on product packaging.

### Inspection of Package Contents

Examine the package carefully before opening to confirm its integrity and the use-by date has not passed. Do not use the device if the package is damaged or opened or if the use-by date has passed.

### Setup

#### Caution

The Emprint™ Ablation Generator with Thermosphere™ Technology, Emprint™ HP Ablation Generator with Thermosphere™ Technology, Emprint™ Ablation Pump, Emprint™ Percutaneous Antennas with Thermosphere™ Technology, and (optional) Emprint™ Ablation Cart, or Emprint™ HP Ablation Cart are designed for use as a single system. Failure to understand and follow instructions provided may result in improper functioning of the system and cause injury to the patient or user.

#### Important

Refer to the system user guide for detailed instructions regarding proper setup and activation of the microwave generator and pump. If using the Emprint™ generator, refer to the *Emprint™ Ablation System with Thermosphere™ Technology User's Guide*. If using the Emprint™ HP generator, refer to the *Emprint™ HP Ablation System with Thermosphere™ Technology User's Guide*.

1. Set up the generator and pump according to the system user's guide.


#### Warning

Do not activate microwave energy until hardware setup is complete and the green, radiating section is fully inserted into tissue.


2. Select the appropriate antenna length (15 cm, 20 cm, or 30 cm) for a given procedure. The shortest length, that enables access to the targeted tissue

for ablation during the planned procedure type, is suggested.

3. Use the appropriate aseptic techniques to remove the sterile antenna, adjustable depth guide, and the cable S clip from the packaging.
4. Open the transparent pump-head cover and insert the orange section of tubing into the pump head. For proper direction of the flow, the arrows on the tubing must point in the same direction as the arrows on the pump head.
5. Close the pump-head cover.
6. Remove the cap from the IV bag spike.
7. Fully insert the bag spike into the port on a one liter or larger sterile saline IV bag. Confirm the IV bag spike extends well into the IV bag's main reservoir.

8. Prime the drip chamber by gently squeezing the sides of the chamber until it is approximately half full of saline. The blue flow indicator should be floating on the surface of the saline.
9. Attach the antenna cable connector to the reusable cable connector. See figure below. The S clip may be attached to the reusable cable to aid in cable management.

- ① Reusable cable connector
- ② Arrows on connectors showing proper alignment between reusable cable and antenna
- ③ Antenna cable connector

### Important

When the arrows are aligned between the antenna cable connect and the reusable cable connector, mate the connectors by pressing them together until an audible click occurs. Rotation of the connectors during connection may cause damage to the device.

### Important

If the generator produces audible alerts and/or illuminates the alert symbols at any time, refer to the troubleshooting instructions in the system user's guide for additional information.

10. Activate the pump by pressing the power switch on the pump's front panel.
11. Allow the pump to fully prime the antenna with fluid before continuing with the procedure. This should take about 30 seconds. The antenna is fully primed when there are no more air bubbles observed circulating through the fluid tubing or the IV bag.

#### Caution

If a leak is encountered, turn OFF the pump and resolve the leak before continuing the procedure.

Verify that the bag spike is securely inserted into the saline bag.

If any other portion of the setup is leaking, replace the antenna as appropriate.

Ensure that the saline flowing to the antenna is constantly dripping into the drip chamber. The pump circulates saline through the antenna and back into the IV bag. Inadequate or no fluid flow can result in unintended thermal injury along the antenna insertion site. If the pump is active and the fluid is not flowing, the tubing may be damaged or blocked.

- Verify the tubing is inserted correctly into the pump
- Verify the arrows on the tubing match the direction of the arrows on the pump
- Verify the bag spike is inserted well into the main reservoir of the IV bag
- Verify the pump has fully primed the antenna and no air bubbles remain in the tubing; if an air bubble appears to be stuck, gently tap at that location to help the trapped air circulate through the system

Retry the procedure. If fluid is still not flowing, replace the antenna.

12. Set the generator's power setting and duration time. (See *Ablation Zone Examples* on page 6 for guidance on power/time settings.)

#### Important

If power and time do not adjust, verify proper connections between the antenna, reusable cable, and generator.

### During a Procedure

1. Remove the protective tube from the antenna shaft.
2. If desired, place the adjustable depth guide onto the antenna shaft. Position the depth guide at the desired depth marker to provide a visual indicator of insertion depth through the skin.

### Warning

The adjustable depth guide is intended for percutaneous use only. This component is not designed to be radio opaque and, therefore, is not intended for use in open procedures.

3. Use image guidance to identify the desired entry location. To prevent damage to the antenna, it is highly recommended that a small skin incision be made prior to insertion. Insert the radiating section of the antenna into the tissue until it is at the intended target site.

#### Warning

Always ensure the microwave energy is inactive before inserting the antenna into the target tissue.

If using an accessory device to aid in the insertion of an antenna (example: cannula or needle guide), ensure:

- The antenna shaft and tip are inserted through the accessory device without resistance or lateral force
- The antenna shaft does not bend or scrape against accessory device
- The accessory device does not apply lateral force to the antenna
- The green, radiating section does not touch bare metal while energy is being delivered

Failure to do so may:

- Cause damage to the antenna including breakage or detachment of the shaft and/or tip
- Cause unintended thermal injury or alter the ablation pattern

The tip of the antenna green, radiating section is sharp. Handle carefully at all times.

Use a scalpel or electrosurgical pencil to expand the insertion point if you meet resistance. Although the antenna is flexible, using force to break through an obstruction may lead to breakage and possible injury to the patient.

Use caution when attempting to ablate lesions close to the surface of the liver. Thermal damage resulting in hemorrhage can occur to adjacent structures.

Do not apply excessive lateral or rotational force during insertion or extraction of the antenna. Doing so may lead to breakage of the antenna and injury to the patient or user.

Do not attempt to bend or reshape the antennas. Product damage and/or malfunction may result.

View the inserted portion of the antenna using image guidance, such as computerized tomography or ultrasound, to ensure proper placement in target tissue and that the antenna is not bent or broken.

### Caution

Confirm proper power and time settings before continuing with the procedure. Always use the shortest time settings that will achieve the desired effect.

If the generator's high-temperature alert triggers during a procedure, turn off the pump, replace the saline bag, and turn on the pump before continuing. After replacing the saline bag, the system resumes at the paused time and power.

### Important

During a procedure, it is important that the antenna remain in the appropriate position throughout the entire procedure. Either manually hold the antenna or affix the antenna in the desired position. If needed, affix the clips provided on the antenna's microwave cable to a stationary object to prevent the cable and antenna from moving during the procedure. Use inter-procedural imaging to monitor placement.

4. Press the generator Start/Stop button or press the footswitch to activate microwave energy to the antenna.

### Important

If a symbol illuminates on the generator, refer to the system user's guide for symbol descriptions and instructions.

Use appropriate image guidance to monitor the formation of the ablation zone.

5. When energy delivery to the antenna stops, allow the pump to continue running for five seconds to cool the antenna. Then rotate the antenna shaft to release any adhered tissue. Then reposition or remove the antenna.

### Warning

Allow the antenna to cool as instructed in step 5 before repositioning or removing the antenna. This reduces the risk of antenna tip detachment.

6. If performing a track ablation, perform this procedure:
  - a) Use the 100 W power output setting whenever possible.
  - b) Set the time value long enough to complete the track ablation.
  - c) Gently rotate the shaft to ensure the antenna can move freely.
  - d) Start microwave energy and begin retracting the antenna from the target tissue at a rate of 0.5 cm per 5 seconds.
  - e) Stop the track ablation when the antenna has been retracted far enough so only the green, radiating section remains in the target tissue.
  - f) Deactivate microwave energy before removing the antenna from the target tissue and patient.

### Warning

Ensure the radiating section is always inserted in tissue when microwave energy is activated. While track ablation is a known technique in RF Ablation, caution should be used in any attempt to cauterize the track using MW ablation. Thermal injury can occur due to radiating field effects.

The green, radiating section may remain hot after activation. Do not touch during or immediately after the application of energy. Do not place the antenna on or near the patient after extraction. Keep all equipment away from flammable materials to avoid combustion.

7. If the antenna is to be reinserted into tissue, clean the green, radiating section by wiping with sterile gauze dampened with sterile saline.

### Note

Do not clean the green, radiating section with any sharp or abrasive objects. Damage to the antenna may occur resulting in diminished performance and possible injury to the patient.

### After a Procedure

1. Verify that the generator and pump are OFF.
2. Open the pump-head cover to equalize pressure throughout the tubing.
3. Disconnect the antenna from the reusable cable and pump.
4. Dispose of the antenna in accordance with your institution's procedures.

### Warning

This product is designed as a single-use device. It has been evaluated in testing to simulate single-procedure conditions. The product is not subjected to reprocessing conditions to evaluate the product's functionality after reprocessing or multiple procedures. Subjecting the product to reprocessing or multiple procedures could potentially impact the device's ability to function as intended.

### Caution

The Emprint™ antenna should be considered a contaminated sharp and it should be disposed of accordingly.

### Ablation Zone Examples

The examples of ablation results on page 83 through page 91 were developed using the Emprint™ Percutaneous Antenna with Thermosphere™ Technology in animal tissue models. Table 1 is from *in-vivo* porcine liver. Table 2 is from *in-vivo* porcine lung. Table 3 is from *in-vivo* porcine kidney. Table 4 is from *ex-vivo* bovine liver. Table 5 is from *ex-vivo* bovine

lung. Table 6 is from *ex-vivo* porcine kidney. The size of ablation zones may vary in clinical settings.

The shaded region in the figure represents the zone of tissue coagulation created around the antenna radiating section.

### More Information

*The Emprint™ Ablation System with Thermosphere™ Technology User's Guide* and *Emprint™ HP Ablation System with Thermosphere™ Technology User's Guide* provide additional information about the Emprint™ Percutaneous Antenna with Thermosphere Technology. Refer to the user's guide for information regarding:

- Using the device
- Other component functions
- Disposal
- Troubleshooting
- Technical specifications
- Technical support

### Antenna Specifications

Temperature Range	Operating: +10°C to +30°C Transport and storage: -30°C to +60°C
Relative Humidity	Operating: 15% to 90% (non-condensing) Transport and storage: 15% to 90% (non-condensing)
Atmospheric Pressure	Operating: 66 kPa to 106 kPa Transport and storage: 50 kPa to 106 kPa