

Shiley™ Flexible Adult Tracheostomy Tubes*

CLINICIAN-INSPIRED SOLUTIONS



*Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.


[Next »](#)



COVIDIEN

positive results for life™

C O N T E N T S

[I. Overview »](#)[II. Features »](#)[III. Ordering Information »](#)

Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.

O V E R V I E W



Design Advances Based on Clinician Input

A focus on care, convenience and comfort

Among other advances, Shiley™ flexible tracheostomy tubes* feature a unique cuff design, color-coding by size (for disposable inner cannula products only), an extended range of sizes, DEHP-Free-plasticizer, clear flange and a design that makes it possible to ventilate with or without an inner cannula. Enhancements like these are focused on three key requirements: patient care, clinical convenience and patient comfort.

*Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.

F E A T U R E S

CARE

Unique taper-shaped cuff

Shiley™ flexible tracheostomy tubes are available with the unique TaperGuard™ taper-shaped cuff.* Compared to the cuffed Shiley™ DCT tracheostomy tube, the Shiley™ flexible tube with TaperGuard™ cuff* can provide a better air and fluid seal at equivalent cuff pressures.¹ And with its low profile, the cuff is designed to reduce insertion and removal force.²

**New material**

Shiley™ flexible tracheostomy tubes were designed to meet the latest safety standards. They are also manufactured with a DEHP-Free-plasticizer.

Enhanced flange

These tubes have a soft recessed flange that is designed to help reduce contact with the stoma site and increase patient comfort. This clear flange allows for enhanced visualization of the underlying skin and stoma site.





Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.

F E A T U R E S

CONVENIENCE

**Percutaneous-capable**

All Shiley™ flexible tracheostomy tubes* have a beveled tip for use with percutaneous-capable loading dilators.

	6.5 mm
	7.0 mm
	7.5 mm
	8.0 mm
	8.5 mm
	9.0 mm
	10.0 mm

Product sizes

Shiley™ flexible tracheostomy tubes* are available in an expanded range of sizes to enable clinicians to find the best and most comfortable fit for different patients.

Color coding (for disposable inner cannula products only)

Shiley™ tracheostomy tubes are now color-coded by size to help clinicians quickly identify the right inner cannula for the corresponding tube.

Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.

F E A T U R E S

COMFORT

Softer material

The new Shiley™ flexible tracheostomy tubes* are made with a softer material than that used in the earlier generation of Shiley tracheostomy tubes.²

Inner cannula

The design of Shiley™ flexible tracheostomy tubes makes it possible to ventilate with or without an inner cannula. Ventilation without the inner cannula effectively increases the diameter of the tube, which may help reduce the work of breathing during the weaning process.³



Shiley™ flexible adult tracheostomy tubes not available in all markets. Shiley™ flexible adult tracheostomy tubes are CE marked and commercially available in European Union countries. Please refer to local product labeling.

ORDERING INFORMATION

**Shiley™ Flexible Adult Tracheostomy Tube with TaperGuard® Cuff***

Disposable Inner Cannula Products		Reusable Inner Cannula Products	Dimensions			
Trach. Tube SKU#	Inner Cannula SKU#	Trach. Tube SKU#	I.D.	O.D.	Length	Inner Cannula I.D.
4CN65A	4IC65	4CN65R	6.5 mm	9.4 mm	62 mm	5.5 mm
5CN70A	5IC70	5CN70R	7.0 mm	10.1 mm	68 mm	6.0 mm
6CN75A	6IC75	6CN75R	7.5 mm	10.8 mm	74 mm	6.5 mm
7CN80A	7IC80	7CN80R	8.0 mm	11.4 mm	77 mm	7.0 mm
8CN85A	8IC85	8CN85R	8.5 mm	12.2 mm	79 mm	7.5 mm
9CN90A	9IC90	9CN90R	9.0 mm	12.7 mm	79 mm	8.0 mm
10CN10A	10IC10	10CN10R	10.0 mm	13.8 mm	79 mm	9.0 mm

Shiley™ Flexible Adult Tracheostomy Tube Cuffless*

Disposable Inner Cannula Products		Reusable Inner Cannula Products	Dimensions			
Trach. Tube SKU#	Inner Cannula SKU#	Trach. Tube SKU#	I.D.	O.D.	Length	Inner Cannula I.D.
4UN65A	4IC65	4UN65R	6.5 mm	9.4 mm	62 mm	5.5 mm
5UN70A	5IC70	5UN70R	7.0 mm	10.1 mm	68 mm	6.0 mm
6UN75A	6IC75	6UN75R	7.5 mm	10.8 mm	74 mm	6.5 mm
7UN80A	7IC80	7UN80R	8.0 mm	11.4 mm	77 mm	7.0 mm
8UN85A	8IC85	8UN85R	8.5 mm	12.2 mm	79 mm	7.5 mm
9UN90A	9IC90	9UN90R	9.0 mm	12.7 mm	79 mm	8.0 mm
10UN10A	10IC10	10UN10R	10.0 mm	13.8 mm	79 mm	9.0 mm



TO CONTACT US, PLEASE VISIT